

ECG SIGNAL FEATURE EXTRACTION

```
IMPORT NUMPY AS NP
IMPORT MATPLOTLIB.PYPLOT AS PLT
FROM SCIPY.FFT IMPORT FFT
FROM SCIPY.SIGNAL IMPORT FIND_PEAKS, BUTTER, FILTFILT

DEF BANDPASS_FILTER(SIGNAL, FS, LOWCUT=0.5, HIGHCUT=50,
ORDER=4):
    """
    APPLIES A BANDPASS FILTER TO REMOVE NOISE FROM THE ECG
    SIGNAL.

    PARAMETERS:
        SIGNAL (ARRAY-LIKE): THE RAW ECG SIGNAL.
        FS (INT): SAMPLING FREQUENCY.
        LOWCUT (FLOAT): LOWER CUTOFF FREQUENCY IN HZ.
        HIGHCUT (FLOAT): HIGHER CUTOFF FREQUENCY IN HZ.
        ORDER (INT): ORDER OF THE BUTTERWORTH FILTER.

    RETURNS:
        ARRAY-LIKE: THE FILTERED SIGNAL.
    """
    NYQUIST = 0.5 * FS
    LOW = LOWCUT / NYQUIST
    HIGH = HIGHCUT / NYQUIST
    B, A = BUTTER(ORDER, [LOW, HIGH], BTYPE='BAND')
    RETURN FILTFILT(B, A, SIGNAL)

DEF EXTRACT_ECG_FEATURES(SIGNAL, FS):
    """
    EXTRACTS KEY FEATURES FROM AN ECG SIGNAL.
```

PARAMETERS:

SIGNAL (ARRAY-LIKE): THE ECG SIGNAL.

FS (INT): SAMPLING FREQUENCY.

RETURNS:

DICT: A DICTIONARY CONTAINING EXTRACTED FEATURES.

"""

FEATURES = {}

TIME-DOMAIN FEATURES

FEATURES['MEAN'] = NP.MEAN(SIGNAL)

FEATURES['STD_DEV'] = NP.STD(SIGNAL)

FEATURES['RMS'] = NP.SQRT(NP.MEAN(NP.SQUARE(SIGNAL)))

FEATURES['ZERO_CROSSINGS'] =

NP.COUNT_NONZERO(NP.DIFF(NP.SIGN(SIGNAL)))

FEATURES['ENERGY'] = NP.SUM(NP.SQUARE(SIGNAL))

DETECT R-PEAKS

**PEAKS, _ = FIND_PEAKS(SIGNAL, HEIGHT=NP.MAX(SIGNAL) * 0.5,
DISTANCE=FS//2)**

FEATURES['NUM_R_PEAKS'] = LEN(PEAKS)

IF LEN(PEAKS) > 1:

RR_INTERVALS = NP.DIFF(PEAKS) / FS

FEATURES['MEAN_RR'] = NP.MEAN(RR_INTERVALS)

FEATURES['STD_RR'] = NP.STD(RR_INTERVALS)

FEATURES['HEART_RATE'] = 60 / NP.MEAN(RR_INTERVALS)

ELSE:

FEATURES['MEAN_RR'] = 0

FEATURES['STD_RR'] = 0

FEATURES['HEART_RATE'] = 0

FREQUENCY-DOMAIN FEATURES USING FFT

N = LEN(SIGNAL)

```

FREQS = NP.FFT.FFTFREQ(N, D=1/FS)
FFT_VALUES = FFT(SIGNAL)
FFT_MAGNITUDE = NP.ABS(FFT_VALUES)

# SPECTRAL FEATURES
SPECTRAL_CENTROID = NP.SUM(FREQS * FFT_MAGNITUDE) /
NP.SUM(FFT_MAGNITUDE)
FEATURES['SPECTRAL_CENTROID'] = SPECTRAL_CENTROID
SPECTRAL_BANDWIDTH = NP.SQRT(NP.SUM(((FREQS -
SPECTRAL_CENTROID)**2) * FFT_MAGNITUDE) /
NP.SUM(FFT_MAGNITUDE))
FEATURES['SPECTRAL_BANDWIDTH'] = SPECTRAL_BANDWIDTH
FFT_MAGNITUDE_NORM = FFT_MAGNITUDE /
NP.SUM(FFT_MAGNITUDE)
SPECTRAL_ENTROPY = -NP.SUM(FFT_MAGNITUDE_NORM *
NP.LOG2(FFT_MAGNITUDE_NORM + 1E-10))
FEATURES['SPECTRAL_ENTROPY'] = SPECTRAL_ENTROPY

RETURN FEATURES

DEF PLOT_ECG_SIGNALS(RAW_SIGNAL, PROCESSED_SIGNAL,
FINAL_SIGNAL, FS):
    """
    PLOTS THREE DIFFERENT ECG SIGNALS: RAW, PROCESSED, AND
    FINAL OUTPUT.

    PARAMETERS:
        RAW_SIGNAL (ARRAY-LIKE): ORIGINAL ECG SIGNAL.
        PROCESSED_SIGNAL (ARRAY-LIKE): FILTERED ECG SIGNAL.
        FINAL_SIGNAL (ARRAY-LIKE): FEATURE-EXTRACTED SIGNAL.
        FS (INT): SAMPLING FREQUENCY.
    """
    N = LEN(RAW_SIGNAL)
    T = NP.ARANGE(0, N) / FS

```

```

PLT.FIGURE(FIGSIZE=(12, 8))

# PLOT RAW SIGNAL
PLT.SUBPLOT(3, 1, 1)
PLT.PLOT(T, RAW_SIGNAL, COLOR="GRAY", LABEL="RAW ECG
SIGNAL")
PLT.TITLE("❑RAW ECG SIGNAL (INPUT)")
PLT.XLABEL("TIME [S]")
PLT.YLABEL("AMPLITUDE")
PLT.LEGEND()

# PLOT PROCESSED SIGNAL
PLT.SUBPLOT(3, 1, 2)
PLT.PLOT(T, PROCESSED_SIGNAL, COLOR="BLUE",
LABEL="FILTERED ECG SIGNAL")
PLT.TITLE("❑PROCESSED ECG SIGNAL (AFTER NOISE
REMOVAL)")
PLT.XLABEL("TIME [S]")
PLT.YLABEL("AMPLITUDE")
PLT.LEGEND()

# PLOT FINAL OUTPUT SIGNAL
PLT.SUBPLOT(3, 1, 3)
PLT.PLOT(T, FINAL_SIGNAL, COLOR="RED", LABEL="FINAL
EXTRACTED FEATURES")
PLT.TITLE("❑FINAL OUTPUT SIGNAL (FEATURE-EXTRACTED)")
PLT.XLABEL("TIME [S]")
PLT.YLABEL("AMPLITUDE")
PLT.LEGEND()

PLT.TIGHT_LAYOUT()
PLT.SHOW()

# EXAMPLE USAGE

```

```

IF __NAME__ == "__MAIN__":
    # GENERATE A SYNTHETIC ECG-LIKE SIGNAL
    FS = 250 # SAMPLING FREQUENCY (ECG SIGNALS TYPICALLY
    250-500 Hz)
    T = NP.Linspace(0, 10, FS * 10) # 10-SECOND ECG SIGNAL
    RAW_SIGNAL = NP.SIN(2 * NP.PI * 1.2 * T) + 0.5 *
    NP.RANDOM.RANDN(LEN(T)) # SIMULATED ECG WITH NOISE

    # APPLY BANDPASS FILTER TO CLEAN THE SIGNAL
    PROCESSED_SIGNAL = BANDPASS_FILTER(RAW_SIGNAL, FS)


    # EXTRACT FEATURES (FINAL SIGNAL BASED ON R-PEAKS)
    FEATURES = EXTRACT_ECG_FEATURES(PROCESSED_SIGNAL, FS)
    FINAL_SIGNAL = NP.ZEROS_LIKE(PROCESSED_SIGNAL)
    PEAKS, _ = FIND_PEAKS(PROCESSED_SIGNAL,
    HEIGHT=NP.MAX(PROCESSED_SIGNAL) * 0.5, DISTANCE=FS//2)
    FINAL_SIGNAL[PEAKS] = PROCESSED_SIGNAL[PEAKS] # SHOW
    ONLY DETECTED R-PEAKS

    # PRINT EXTRACTED FEATURES
    PRINT("EXTRACTED ECG FEATURES:")
    FOR KEY, VALUE IN FEATURES.ITEMS():
        PRINT(F"{KEY}: {VALUE:.5F}")


    # PLOT ALL THREE SIGNALS
    PLOT_ECG_SIGNALS(RAW_SIGNAL, PROCESSED_SIGNAL,
    FINAL_SIGNAL, FS)

```


1 Raw ECG Signal (Input)



2 Processed ECG Signal (After Noise Removal)



3 Final Output Signal (Feature-Extracted)



ECG SIGNAL ABNORMALITY DETECTION

```
import numpy as np
import matplotlib.pyplot as plt
from scipy.signal import find_peaks, butter, filtfilt

def bandpass_filter(signal, fs, lowcut=0.5, highcut=50,
                    order=4):
    """
    Applies a bandpass filter to remove noise from the ECG
    signal.
    """
    nyquist = 0.5 * fs
    low = lowcut / nyquist
    high = highcut / nyquist
    b, a = butter(order, [low, high], btype='band')
    return filtfilt(b, a, signal)

def detect_abnormalities(signal, fs):
    """
    Detects abnormalities in the ECG signal based on heart
    rate variations.
    """
    peaks, _ = find_peaks(signal,
                          height=np.percentile(signal, 95), distance=fs//2)

    if len(peaks) < 2:
        return np.zeros_like(signal), 0, 0, "No R-peaks
        detected"

    rr_intervals = np.diff(peaks) / fs
```

```
AVG_RR = NP.MEAN(RR_INTERVALS)
STD_RR = NP.STD(RR_INTERVALS)
HEART_RATE = 60 / AVG_RR IF AVG_RR > 0 ELSE 0
```

```
ABNORMAL_SIGNAL = NP.ZEROS_LIKE(SIGNAL)
FOR I, RR IN ENUMERATE(RR_INTERVALS):
    IF ABS(RR - AVG_RR) > 1.5 * STD_RR:
        ABNORMAL_SIGNAL[PEAKS[I]] = SIGNAL[PEAKS[I]]
```

```
# CLASSIFY ABNORMALITY
IF HEART_RATE > 100:
    CONDITION = "TACHYCARDIA (FAST HEART RATE)"
ELIF HEART_RATE < 60:
    CONDITION = "BRADYCARDIA (SLOW HEART RATE)"
ELSE:
    CONDITION = "NORMAL ECG"
```

```
RETURN ABNORMAL_SIGNAL, HEART_RATE, STD_RR, CONDITION
```

```
DEF PLOT_ECG_SIGNALS(RAW_SIGNAL, PROCESSED_SIGNAL,
ABNORMAL_SIGNAL, FS):
```

```
    """
```

```
    PLOTS THREE ECG SIGNALS: RAW, PROCESSED, AND
    ABNORMALITY DETECTION.
```

```
    """
```

```
    N = LEN(RAW_SIGNAL)
```

```
    T = NP.ARANGE(N) / FS
```

```
    PLT.FIGURE(FIGSIZE=(12, 8))
```

```
    PLT.SUBPLOT(3, 1, 1)
```

```
    PLT.PLOT(T, RAW_SIGNAL, COLOR="GRAY", LABEL="RAW ECG
    SIGNAL")
```

```
    PLT.TITLE("RAW ECG SIGNAL")
```

```
    PLT.XLABEL("TIME [S]")
```



```
PLT.YLABEL("AMPLITUDE")
PLT.LEGEND()
```

```
PLT.SUBPLOT(3, 1, 2)
PLT.PLOT(T, PROCESSED_SIGNAL, COLOR="BLUE",
LABEL="FILTERED ECG SIGNAL")
PLT.TITLE("PROCESSED ECG SIGNAL")
PLT.XLABEL("TIME [S]")
PLT.YLABEL("AMPLITUDE")
PLT.LEGEND()
```

```
PLT.SUBPLOT(3, 1, 3)
PLT.PLOT(T, PROCESSED_SIGNAL, COLOR="BLUE", ALPHA=0.5,
LABEL="FILTERED ECG")
PLT.SCATTER(T[ABNORMAL_SIGNAL > 0],
ABNORMAL_SIGNAL[ABNORMAL_SIGNAL > 0], COLOR="RED",
LABEL="ABNORMALITIES", ZORDER=3)
PLT.TITLE("DETECTED ABNORMALITIES (RED DOTS)")
PLT.XLABEL("TIME [S]")
PLT.YLABEL("AMPLITUDE")
PLT.LEGEND()
```

```
PLT.TIGHT_LAYOUT()
PLT.SHOW()
```

EXAMPLE USAGE


```
IF __NAME__ == "__MAIN__":
    FS = 250 # SAMPLING FREQUENCY
    T = NP.Linspace(0, 10, FS * 10) # 10-SECOND ECG SIGNAL
    RAW_SIGNAL = NP.SIN(2 * NP.PI * 1.2 * T) + 0.5 *
NP.RANDOM.RANDN(LEN(T)) # SIMULATED ECG WITH NOISE

    PROCESSED_SIGNAL = BANDPASS_FILTER(RAW_SIGNAL, FS)
    ABNORMAL_SIGNAL, HEART_RATE, STD_RR, CONDITION =
DETECT_ABNORMALITIES(PROCESSED_SIGNAL, FS)
```


```
PRINT(F"HEART RATE: {HEART_RATE:.2F} BPM")  
PRINT(F"RR INTERVAL VARIABILITY: {STD_RR:.2F}")  
PRINT(F"CONDITION: {CONDITION}")
```

```
PLOT_ECG_SIGNALS(RAW_SIGNAL, PROCESSED_SIGNAL,  
ABNORMAL_SIGNAL, FS)
```

Raw ECG Signal



Processed ECG Signal



Detected Abnormalities (Red Dots)

