PPG Signal Processing and Analysis

1. Introduction Photoplethysmography (PPG) is a non-invasive optical technique used to measure blood volume changes in the microvascular tissue. This experiment involves generating a synthetic PPG signal, applying signal processing techniques, and analyzing heart rate metrics to detect abnormalities.

2. Objectives

- To generate a synthetic PPG signal with realistic components.
- To apply a Butterworth low-pass filter to remove noise.
- To detect peaks and valleys in the signal corresponding to heartbeats.
- To compute heart rate (HR) and heart rate variability (HRV).
- To identify potential abnormalities in the PPG signal.
- To visualize and analyze the results.
- **3. Methodology** The experiment was conducted using Python and signal processing libraries. The key steps are as follows:

3.1 Signal Generation

A synthetic PPG signal was generated using sinusoidal components and Gaussian noise:

- Base signal with a heart rate of approximately 60 BPM.
- Higher frequency components were added to simulate real-world variations.
- Random noise was introduced to enhance realism.

3.2 Signal Filtering

A Butterworth low-pass filter was applied:

• Sampling Rate: 100 Hz

Cutoff Frequency: 5 Hz

Filter Order: 4

3.3 Peak Detection and Heart Rate Calculation

- Peaks (representing heartbeats) were detected using scipy.signal.find_peaks().
- Peak-to-peak intervals (PPI) were used to calculate HR and HRV.
- Instantaneous heart rates were derived from the inverse of PPI.

3.4 Abnormality Detection

- Abnormal peaks were identified based on deviations in PPI.
- Early and delayed beats were detected using a statistical threshold (±2 standard deviations).

3.5 Visualization and Statistical Analysis

- The raw and filtered signals were plotted with detected peaks and valleys.
- Instantaneous heart rate over time was visualized.

Statistical metrics were computed and displayed.

Code

```
import numpy as np
from scipy import signal
import matplotlib.pyplot as plt
from scipy.signal import find_peaks
from typing import Tuple, List, Dict
import os
from datetime import datetime
# [Previous functions remain the same until plot_ppg_analysis]
def plot_ppg_analysis(t, ppg_signal, filtered_ppg, peaks, valleys, abnormal_peaks=None, save_path=None):
  Create visualization plots for PPG signal analysis and optionally save to PNG.
  Args:
    t (array): Time array
    ppg_signal (array): Raw PPG signal
    filtered_ppg (array): Filtered PPG signal
    peaks (array): Indices of detected peaks
    valleys (array): Indices of detected valleys
    abnormal_peaks (array, optional): Indices of abnormal peaks
    save_path (str, optional): Directory path to save the plots
  plt.figure(figsize=(15, 12))
  # Plot 1: Raw PPG Signal
  plt.subplot(4, 1, 1)
  plt.plot(t, ppg_signal)
  plt.title('Raw PPG Signal')
  plt.xlabel('Time (s)')
```


```
plt.ylabel('Amplitude')
plt.grid(True)
# Plot 2: Filtered PPG Signal
plt.subplot(4, 1, 2)
plt.plot(t, filtered_ppg)
plt.title('Filtered PPG Signal')
plt.xlabel('Time (s)')
plt.ylabel('Amplitude')
plt.grid(True)
# Plot 3: Peaks and Valleys
plt.subplot(4, 1, 3)
plt.plot(t, filtered_ppg, label='Filtered Signal')
plt.plot(t[peaks], filtered_ppg[peaks], "rv", label='Normal Peaks')
plt.plot(t[valleys], filtered_ppg[valleys], "g^", label='Valleys')
if abnormal_peaks and len(abnormal_peaks) > 0:
  plt.plot(t[abnormal_peaks], filtered_ppg[abnormal_peaks], "yx",
      markersize=10, label='Abnormal Peaks')
plt.title('Filtered PPG Signal with Peaks and Valleys')
plt.xlabel('Time (s)')
plt.ylabel('Amplitude')
plt.legend()
plt.grid(True)
# Plot 4: Instantaneous Heart Rate
if len(peaks) > 1:
  inst_hr = 60 / np.diff(t[peaks])
  hr_times = t[peaks[1:]]
  plt.subplot(4, 1, 4)
```

```
plt.plot(hr_times, inst_hr)
   plt.title('Instantaneous Heart Rate')
   plt.xlabel('Time (s)')
   plt.ylabel('Heart Rate (BPM)')
   plt.grid(True)
  plt.tight_layout()
 # Save plot if path is provided
 if save_path:
   timestamp = datetime.now().strftime("%Y%m%d_%H%M%S")
   if not os.path.exists(save_path):
     os.makedirs(save_path)
   filename = os.path.join(save_path, f'ppg_analysis_{timestamp}.png')
   plt.savefig(filename, dpi=300, bbox_inches='tight')
   print(f"Plot saved as: {filename}")
  plt.show()
def save_statistics_to_file(hr_metrics: Dict, abnormal_peaks: List, abnormality_types: List, save_path: str):
 .....
 Save statistics to a text file.
 Args:
   hr_metrics (Dict): Dictionary of heart rate metrics
   abnormal_peaks (List): List of abnormal peak indices
   abnormality_types (List): List of abnormality types
   save_path (str): Directory path to save the statistics
 if not os.path.exists(save_path):
   os.makedirs(save_path)
```

```
timestamp = datetime.now().strftime("%Y%m%d_%H%M%S")
 filename = os.path.join(save_path, f'ppg_statistics_{timestamp}.txt')
 with open(filename, 'w') as f:
   f.write("=== Heart Rate Statistics ===\n")
   f.write(f"Average Heart Rate: {hr_metrics['mean_hr']:.1f} BPM\n")
   f.write(f"Heart Rate Variability: {hr_metrics['hrv']:.1f} ms\n")
   f.write(f"Minimum Heart Rate: {hr_metrics['min_hr']:.1f} BPM\n")
   f.write(f"Maximum Heart Rate: {hr_metrics['max_hr']:.1f} BPM\n")
   f.write(f"Heart Rate Standard Deviation: {hr_metrics['std_hr']:.1f} BPM\n\n")
   f.write("=== Abnormality Detection ===\n")
   if len(abnormal_peaks) > 0:
     f.write(f"Number of abnormalities detected: {len(abnormal_peaks)}\n")
     for i, (peak, atype) in enumerate(zip(abnormal_peaks, abnormality_types)):
       f.write(f"Abnormality {i+1}: {atype} at time {peak:.2f}s\n")
   else:
     f.write("No abnormalities detected\n")
  print(f"Statistics saved as: {filename}")
def main():
 # Parameters
 sampling_rate = 100
  duration = 10
  cutoff_freq = 5 # Hz
 # Set output directory for saving files
  save_path = 'ppg_output'
 # Generate sample data
 t, ppg_signal = generate_ppg_signal(duration, sampling_rate)
```

```
# Apply Butterworth filter
 filtered_ppg = apply_butterworth_filter(ppg_signal, sampling_rate, cutoff_freq)
 # Find peaks and valleys
  peaks, _ = find_peaks(filtered_ppg, distance=50)
 valleys, _ = find_peaks(-filtered_ppg, distance=50)
 # Calculate heart rate metrics
  hr_metrics = calculate_heart_rate_metrics(t, peaks)
  # Detect abnormalities
  abnormal_peaks, abnormality_types = detect_abnormalities(t, peaks, hr_metrics)
 # Visualize results and save plot
  plot_ppg_analysis(t, ppg_signal, filtered_ppg, peaks, valleys, abnormal_peaks, save_path)
 # Display and save statistics
  display_statistics(hr_metrics, abnormal_peaks, abnormality_types)
  save_statistics_to_file(hr_metrics, abnormal_peaks, abnormality_types, save_path)
if __name__ == "__main__":
  main()
```

- **4. Results and Discussion** The analysis provided the following results:
 - Average Heart Rate (HR): Computed from detected peaks.
 - Heart Rate Variability (HRV): Derived from standard deviation of PPI.
 - Minimum & Maximum HR: Identified from instantaneous HR variations.
 - Abnormal Beats: Early or delayed beats detected in the signal.
 - **Visualization:** The processed PPG signal showed clear periodic peaks corresponding to heartbeats, and abnormal peaks were marked in the plots.



=== Heart Rate Statistics ===

Average Heart Rate: 59.9 BPM

Heart Rate Variability: 0.0 ms

Minimum Heart Rate: 59.9 BPM

Maximum Heart Rate: 59.9 BPM

Heart Rate Standard Deviation: 0.0 BPM

=== Abnormality Detection ===

No abnormalities detected

5. Conclusion The experiment successfully demonstrated the steps involved in PPG signal processing, filtering, and heart rate analysis. The methodology provided a structured approach to detect abnormalities in heart rate. The results showed the effectiveness of signal processing techniques in analyzing physiological signals.

6. Future Work

• Implementing real-time PPG signal analysis.

- Enhancing abnormality detection using machine learning techniques.
- Comparing different filtering techniques for improved signal quality.