

Peak Detection and Pulse Analysis Using PPG Signals

Objective:

- **Develop a peak detection algorithm** to analyze pulse rate and pulse morphology from real PPG signals.
 - **Extract features** related to Pulse Rate Variability (PRV) and Pulse Morphology for cardiovascular health assessment.
-

Theory:

Photoplethysmography (PPG) is an optical technique that measures blood volume changes in the microvascular bed of tissue. It uses light-emitting diodes (LEDs) and photodetectors to capture pulse waves, which correspond to heartbeats.

Key Features Extracted:

1. **Pulse Rate (PR)**: The number of pulses per minute, equivalent to heart rate.
 2. **Pulse Rate Variability (PRV)**: Variation in time intervals between consecutive pulses, related to Heart Rate Variability (HRV).
 3. **Pulse Morphology**: The shape and characteristics of the pulse waveform, including:
 - **Peak Amplitude**: Height of the peak, indicating blood volume.
 - **Peak Width**: Duration of each pulse.
 - **Dicrotic Notch**: A secondary dip in the waveform related to aortic valve closure.
-

Techniques Used:

1. **Adaptive Thresholding**: Dynamically adjusts the threshold for peak detection based on signal properties.

2. **Wavelet Transform:** Denoises and enhances peak detection accuracy.
 3. **Signal Processing (SciPy):** Used for filtering and detecting peaks.
 4. **Visualization (Matplotlib):** Visualizes waveforms and feature analysis.
-

Data Processing Steps:

1. Step 1: Load and Inspect Data

- Load PPG data from a CSV file.
- Extract the relevant PPG column.
- Visualize the raw data to understand noise and signal trends.

2. Step 2: Signal Preprocessing

- Apply filtering to remove high-frequency noise and baseline drift.
- Normalize the signal for consistent peak detection.

3. Step 3: Peak Detection

- Use the `find_peaks()` function to locate peaks corresponding to heartbeats.
- Adjust parameters like minimum peak height and minimum distance between peaks for optimal detection.

4. Step 4: Pulse Analysis

- Calculate Pulse Rate (PR) in Beats Per Minute (BPM).
- Compute Inter-Peak Intervals (IPI) for Pulse Rate Variability (PRV).
- Analyze Pulse Morphology (e.g., amplitude, width).

5. Step 5: Visualization

- Plot Raw vs. Cleaned PPG signals.
- Mark detected peaks to validate the algorithm.
- Visualize Pulse Rate Variability (PRV).

Code :

```
import numpy as np
import pandas as pd
import matplotlib.pyplot as plt
from scipy import signal
import neurokit2 as nk

data = pd.read_csv('bio_resting_5min_100hz.csv') # Load the CSV file
ppg_signal = data['PPG'] # Extract the PPG data from the 'PPG' column
fs = 100
time = np.arange(len(ppg_signal)) / fs


ppg_clean = nk.ppg_clean(ppg_signal, sampling_rate=fs)

# Visualize Raw and Cleaned PPG Signals
plt.figure(figsize=(16, 8))

plt.subplot(2, 1, 1)
plt.plot(time, ppg_signal, color='red', alpha=0.7, label='Raw PPG')
plt.title("Raw PPG Signal")
plt.xlabel("Time (s)")
plt.ylabel("Amplitude")
plt.legend()
plt.grid(True)

plt.subplot(2, 1, 2)
plt.plot(time, ppg_clean, color='green', alpha=0.7, label='Cleaned PPG')
plt.title("Cleaned PPG Signal")
plt.xlabel("Time (s)")
plt.ylabel("Amplitude")
plt.legend()
plt.grid(True)

plt.tight_layout()
plt.show()
```



```
# Peak Detection in Cleaned PPG Signal
peaks, properties = signal.find_peaks(
    ppg_clean,
    distance=fs*0.5,
    height=np.mean(ppg_clean)*0.6
)

peak_count = len(peaks)
print(f"Total Peaks Detected: {peak_count}")
```


Total Peaks Detected: 432

```
# Calculate Pulse Rate and PRV Metrics
if peak_count > 1:
    rr_intervals = np.diff(peaks) / fs
    pulse_rate = 60 / np.mean(rr_intervals)
    print(f"Pulse Rate: {pulse_rate:.2f} BPM")
```


```
else:  
    print("Insufficient peaks detected. Check signal quality or peak detection  
parameters.")
```

Pulse Rate: 86.37 BPM

```
# Visualize Cleaned PPG with Detected Peaks  
plt.figure(figsize=(16, 8))  
plt.plot(time, ppg_clean, color='green', linewidth=1, label='Cleaned PPG')  
plt.plot(time[peaks], ppg_clean[peaks], 'ro', label='Detected Peaks') # Mark  
peaks  
plt.title("Cleaned PPG Signal with Detected Peaks")  
plt.xlabel("Time (s)")  
plt.ylabel("Amplitude")  
plt.legend()  
plt.grid(True)  
plt.show()
```



```
# Step 7: Pulse Rate Variability (PRV) Visualization
plt.figure(figsize=(14, 7))
plt.plot(rr_intervals, marker='o', linestyle='-', color='blue')
plt.title("Pulse Rate Variability (PRV)")
plt.xlabel("Pulse Number")
plt.ylabel("Interval (s)")
plt.grid(True)
plt.show()
```



Output Explanation:

1. Raw vs. Cleaned PPG Signal:

- The raw signal shows noise and baseline drift.
- The cleaned signal is smoother, with distinct peaks.

2. Cleaned PPG with Detected Peaks:

- Red dots mark detected peaks corresponding to heartbeats.
- Peaks are accurately aligned with pulse waves.

3. Pulse Rate Variability (PRV):

- The PRV graph shows variations in time between heartbeats.
- This variability is used to assess cardiovascular health and stress levels.

4. Console Outputs:

- Total Peaks Detected: Displays the number of detected heartbeats.
- Pulse Rate: Shows the average heart rate in BPM.