

Department of Information & Communication Engineering

Course Title: Signals and Systems

Course Code: ICE-2204

Lab Name: Extracting

abnormality from raw PPG signal.

Submitted To:

Dr. Md. Imran Hossain Associate Professor Department of ICE Pabna University of Science &Technology

Submitted By:

Mst. Tahomina Kabir Rany

Roll: 210642 Session: 2021-22 2nd Year 2nd Semester

Dept. of Information & Communication

Engineering

Pabna University of science & Technology

Date of Submission: 14.02.2025

Introduction to PPG (Photoplethysmography)

Photoplethysmography (PPG) is an optical technique used to measure blood volume changes in the microvascular tissue. It is widely employed in wearable health devices and medical monitoring systems to track heart rate (HR), blood oxygen levels (SpO₂), and other cardiovascular parameters.

PPG signals are acquired using light-emitting diodes (LEDs) and photodetectors, which detect changes in light absorption due to pulsatile blood flow.

Extracting Raw PPG Signals

The process of obtaining a raw PPG signal involves several key steps:

Hardware Components

- **Light Source (LEDs):** Typically, red, green, or infrared (IR) LEDs are used. Green light is often preferred for wrist-based devices due to better penetration in blood-rich tissues.
- Photodetector: Measures the intensity of reflected or transmitted light, capturing pulsatile variations in blood volume.

Signal Acquisition

- The raw signal is an optical intensity waveform representing blood volume changes.
- Typically, a sample rate of 100–1000 Hz is used for high-resolution data.

Preprocessing of PPG Signals

Since raw PPG signals contain noise and artifacts, preprocessing is crucial:

Signal Filtering:

- Low-pass filtering (e.g., Butterworth, moving average) to remove high-frequency noise.
- o High-pass filtering to eliminate baseline drift.
- Notch filtering to remove power line interference (50/60 Hz).

Motion Artifact Removal:

- Adaptive filtering or signal decomposition techniques like wavelet transform, Empirical Mode Decomposition (EMD).
- Accelerometer-based correction for wearable devices.

Normalization & Smoothing:

- Normalization to scale the signal amplitude.
- Smoothing using Savitzky-Golay filters or median filtering.

Detecting Abnormalities in PPG Signals

PPG signals can be analyzed to detect cardiovascular abnormalities using various signal processing and machine learning techniques.

Heart Rate and Heart Rate Variability (HRV) Analysis

- Extracting peaks from the PPG waveform provides heart rate (HR).
- HRV analysis (time and frequency domain) helps detect arrhythmias and autonomic dysfunction.

Abnormality Detection Techniques

1. Arrhythmia Detection:

- Irregular peak intervals indicate atrial fibrillation (AF) or bradycardia.
- Machine learning classifiers (e.g., Random Forest, CNNs) can detect arrhythmias from PPG waveforms.

2. Blood Pressure Estimation:

 Pulse Transit Time (PTT) and machine learning models estimate systolic and diastolic BP.

3. SpO₂ (Oxygen Saturation) Estimation:

- Dual-wavelength PPG (red & IR) is used.
- Abnormally low SpO₂ (<90%) may indicate hypoxia or respiratory distress.

4. Respiratory Rate (RR) Estimation:

 Respiratory-induced modulations in PPG amplitude and frequency are analyzed.

5. **Detecting Shock or Hypovolemia:**

 Decreased amplitude and variations in PPG signal morphology may indicate shock or hypovolemia.

Machine Learning and AI in PPG Abnormality Detection

With advancements in AI, deep learning-based approaches improve accuracy:

- Convolutional Neural Networks (CNNs): Extract features from PPG signals.
- **Recurrent Neural Networks (RNNs, LSTMs):** Capture temporal dependencies in time-series data.

Source Code: Given below

```
import numpy as np
import scipy.signal as signal
import matplotlib.pyplot as plt
from scipy.signal import find_peaks, butter, filtfilt
```


```
# Generate synthetic PPG signal (replace this with your raw PPG data)
def generate_synthetic_ppg(length=1000, sampling_rate=100):
    t = np.linspace(0, length / sampling_rate, length)
    ppg = 0.5 * np.sin(2 * np.pi * 1.2 * t) + 0.05 *
np.random.randn(length)
    return t, ppg
```

```
# Bandpass filter
def bandpass filter(data, lowcut, highcut, fs, order=4):
  nyquist = 0.5 * fs
  low = lowcut / nyquist
  high = highcut / nyquist
  b, a = butter(order, [low, high], btype="band")
  filtered data = filtfilt(b, a, data)
  return filtered data
# Detect peaks (heartbeats)
def detect peaks(ppg, sampling rate, distance=50):
  peaks, = find peaks(ppg, distance=distance)
  return peaks
# Calculate heart rate variability (HRV)
def calculate_hrv(peaks, sampling_rate):
  rr intervals = np.diff(peaks) / sampling rate # Time between
consecutive peaks
  mean rr = np.mean(rr intervals)
  std rr = np.std(rr intervals)
  return rr_intervals, mean_rr, std_rr
```


```
# Abnormality detection based on HRV thresholds
def detect_abnormalities(rr_intervals, mean_rr, std_rr, threshold=0.2):
  abnormalities = []
  for i, rr in enumerate(rr intervals):
    if abs(rr - mean rr) > threshold * std rr:
      abnormalities.append(i)
  return abnormalities
# Main function
def main():
  # Parameters
  sampling rate = 100 # Hz
  lowcut = 0.5
                  # Hz
  highcut = 5.0
                  # Hz
  threshold = 0.2 # Abnormality detection threshold
  # Generate or load PPG data
  t, raw_ppg = generate_synthetic_ppg(length=3000,
sampling rate=sampling rate)
  # Preprocess PPG data
```

```
filtered ppg = bandpass filter(raw ppg, lowcut, highcut,
sampling rate)
  # Detect peaks
  peaks = detect peaks(filtered ppg, sampling rate, distance=int(0.6 *
sampling rate)) # Min distance ~600ms
  # Calculate HRV
  rr intervals, mean rr, std rr = calculate hrv(peaks, sampling rate)
  # Detect abnormalities
  abnormalities = detect abnormalities(rr intervals, mean rr, std rr,
threshold=threshold)
  # Plot results
  plt.figure(figsize=(12, 6))
  plt.plot(t, filtered ppg, label="Filtered PPG")
  plt.scatter(t[peaks], filtered ppg[peaks], color="red", label="Detected
Peaks")
  plt.title("PPG Signal with Detected Peaks")
  plt.xlabel("Time (s)")
  plt.ylabel("Amplitude")
  plt.legend()
```

```
plt.grid()
  plt.figure(figsize=(12, 6))
  plt.plot(rr intervals, label="RR Intervals (s)")
  plt.scatter(abnormalities, np.array(rr intervals)[abnormalities],
color="red", label="Abnormal Intervals")
  plt.title("RR Intervals with Abnormalities")
  plt.xlabel("Interval Index")
  plt.ylabel("RR Interval (s)")
  plt.legend()
  plt.grid()
  plt.show()
  # Print summary
  print(f"Mean RR Interval: {mean rr:.2f} s")
  print(f"Standard Deviation of RR Intervals: {std rr:.2f} s")
  print(f"Number of Abnormalities Detected: {len(abnormalities)}")
if __name__ == "__main__":
  main()
```



☆ ← → | 中 Q ∓ | □ (x, y) = (5.55, 0.586)



☆ ♦ ♦ 4 Q ≨ 🖺