Purpose:

```
#creating a ppg signal with noise.
#Filltered the ppg signal by low pass fillter.
#Detected peaks and valleys.
#Abnormal peaks detection.
# And Count number of abnormal peaks,
estimated heart rate (BPM) and
total peak .
```

code:

```
import numpy as np
import pandas as pd
import matplotlib.pyplot as plt
from scipy.signal import find peaks, butter, filtfilt
# Generate sample PPG signal
np.random.seed(0)
time = np.linspace(0, 10, 1000)
ppq signal = np.sin(2 * np.pi * 1.2 * time) + 0.5 *
np.random.normal(size=len(time))
# Define the lowpass filter function
def butter lowpass filter(data, cutoff, fs, order=5):
    nyquist = 0.5 * fs
    normal cutoff = cutoff / nyquist
    b, a = butter(order, normal cutoff, btype='low',
analog=False)
    y = filtfilt(b, a, data)
    return y
```


```
# Filter settings
fs = 100 # Sampling frequency in Hz
cutoff = 3 # Cutoff frequency in Hz
filtered ppg = butter lowpass filter(ppg signal, cutoff, fs)
# Peak and valley detection
peaks, = find peaks(filtered ppg, height=0.5,
distance=fs//2)
valleys, = find peaks(-filtered ppg, height=0.5,
distance=fs//2)
# Abnormal peaks detection
peak heights = filtered ppg[peaks]
abnormal peaks = peaks[peak heights > 1] # Threshold for
high spikes
# Calculate Heart Rate (BPM)
time diffs = np.diff(time[peaks]) # Time difference between
consecutive peaks
avg rr interval = np.mean(time diffs) # Average R-R interval
in seconds
heart rate = 60 / avg rr interval # Convert to BPM
print(f"Number of abnormal peaks: {len(abnormal peaks)}")
print(f"Estimated Heart Rate: {heart rate:.2f} BPM")
# Count total peak values
total peaks = len(peaks)
print(f"Total Peak Count: {total peaks}")
# Plot 1: Raw PPG Signal
plt.figure(figsize=(12, 6))
```


```
plt.plot(time, ppg signal, label="Raw PPG Signal", alpha=0.5)
plt.title("Raw PPG Signal")
plt.xlabel("Time (s)")
plt.ylabel("Amplitude")
plt.grid()
plt.legend()
plt.show()
# Plot 2: Filtered PPG Signal
plt.figure(figsize=(12, 6))
plt.plot(time, filtered ppg, label="Filtered PPG Signal",
linewidth=2)
plt.title("Filtered PPG Signal")
plt.xlabel("Time (s)")
plt.ylabel("Amplitude")
plt.grid()
plt.legend()
plt.show()
# Plot 3: Peaks and Valleys
plt.figure(figsize=(12, 6))
plt.plot(time, filtered ppg, label="Filtered PPG Signal",
linewidth=2)
plt.plot(time[peaks], filtered ppg[peaks], "go",
label="Detected Peaks")
plt.plot(time[valleys], filtered ppg[valleys], "ro",
label="Detected Valleys")
plt.title("Detected Peaks and Valleys")
plt.xlabel("Time (s)")
plt.ylabel("Amplitude")
plt.grid()
plt.legend()
plt.show()
```


Plot 4: Abnormal Peaks


```
plt.figure(figsize=(12, 6))
plt.plot(time, filtered_ppg, label="Filtered PPG Signal",
  linewidth=2)
plt.plot(time[abnormal_peaks], filtered_ppg[abnormal_peaks],
  "kx", label="Abnormal Peaks", markersize=10)
plt.title("Abnormal Peaks Detection")
plt.xlabel("Time (s)")
plt.ylabel("Amplitude")
plt.grid()
plt.legend()
plt.show()
```

Output:









Number of abnormal peaks: 4

Estimated Heart Rate: 71.98 BPM

Total Peak Count: 12